

PLASTIPAD®

REUSABLE HYPER-HYPOTHEMIA BLANKETS

The PlastiPad® is a reusable water blanket that is durable and convenient. It may resist punctures and has a seam strength that prolongs the life of the blankets. These blankets have a non-porous surface that is easy to clean off blood, soil, and/or debris.

The PlastiPad® blanket utilizes a random flow pattern for uniform distribution of water, yielding fast and efficient circulation throughout the entire blanket to provide effective patient therapy.



PlastiPad® Product Description

Model 196	Adult - 24" x 60" (61cm x 152.4cm) (Not available in Europe)
Model 196CPC	Adult - 24" x 60" (61cm x 152.4cm) (requires #286 connecting hose)
Model 195N	Narrow Adult - 20" x 60" (50.8cm x 152.4cm) (requires #286 connecting hose)
Model 194	Pediatric - 22" x 30" (55.9cm x 76.2cm) (Not available in Europe)
Model 194CPC	Pediatric - 22" x 30" (55.9cm x 76.2cm) (requires #286 connecting hose)
Model 193	Infant - 12" x 18" (30.5cm x 45.7cm) (Not available in Europe)
Model 193CPC	Infant - 12" x 18" (30.5cm x 45.7cm) (requires #286 connecting hose)