

ELECTRI-COOL® II

LOCALIZED COLD THERAPY SYSTEM

The Electri-Cool® II is the premier choice when it comes to localized cold therapy systems. The Electri-Cool® II may be attached to the patient's bed with a bracket to save valuable floor space or may be mounted on a stand for easy mobility. It was designed for ease-of-use, just simply plug in the unit, turn on the power and select your desired level of therapy. The Electri-Cool® II saves valuable time with automatic low water notification when the unit needs to be refilled with water.



Product Description

Model 767	110/115v & 220/240v
Dimensions	12" W x 6" D x 12" H (30.5cm x 15.2cm x 30.5cm)
Water Temp. Range	5°C - 7°C (40°F - 45°F) 8°C - 10°C (45°F - 50°F) 11°C - 13°C (50°F - 55°F) COOLING ONLY