

Normothermia Systems



CSZ's normothermia systems help keep your patient Normothermic before, during and after surgery.

CSZ offers three different systems for providing normothermia; the Norm- O-Temp® is water based, the WarmAir® is an air based system and the SurfaceTemp® which is resistive warming. All of these systems provide the desired warmth to the patient; it just depends on the caregiver's preference.

The Norm-O-Temp® combined with the Gelli-Roll® offers excellent temperature control while providing complete access to the patient. It is a great solution for the Operating Room and Emergency Department.

The WarmAir® has 3 therapeutic temperature settings that help keep the patient at a normothermic level. Any department in the hospital that is looking for an effective active warming system, the WarmAir® is your solution.

The SurfaceTemp® is a resistant warming technology OR Table Pad. It redistributes pressure and provides comfortable consistent and active warming to every patient.

Using these 3 systems together for combination therapy in specific cases can meet the AORN guidance for warming.

