

STAFF VEST

HYPER-HYPOTHEMRIA VEST FOR CLINICAL STAFF

The Staff Cooling Vest helps keep the surgical team comfortable throughout a procedure and is designed to circulate water evenly and effectively without the mess of ice. With the Blanketrol® II or III, surgeons and their staff can pre-condition and control the water temperature for warming or cooling.

- Sterile gowns and/or lead vest may be worn on top of the Staff Vest
- One size fits most staff members
- Adjustable Velcro® Straps
- Can be connected to CSZ equipment below the sterile field
- Lightweight and permanently attached hoses



Staff Vest Product Description

Model 300

one size fits most