


FILTEREDFLO® BLANKETS

CONVECTIVE WARMING BLANKETS

CSZ uses Filtered Air Warming as opposed to Forced Air Warming.

The unique patented design of FilteredFlo® blankets permits use of a lower velocity blower to supply gently moving, filtered air. The blankets also act as a second filter to the air moving through the blower before it reaches the patient.



Product Description

Model 243	Adult Full Body Blanket 92" x 50" (233.7cm x 127cm)
Model 462	Adult Full Body Blanket 80" x 36" (203cm x 91.4cm)
Model 443	Upper Body Blanket 88" x 30" (223.5cm x 76.2cm)
Model 461	Upper Body Blanket 76" x 24" (193cm x 60.9cm)
Model 442	Lower Body Blanket 53" x 40" (134.6cm x 101.6cm)
Model 460	Lower Body Blanket 46" x 36" (116.8cm x 91.4cm)
Model 244	Pediatric Blanket 56" x 40" (142.2cm x 101.6cm)
Model 246	Infant Blanket (Under/Over) 24" x 37" (61cm x 94cm)
Model 248	Adult Underbody Blanket 78.5" x 40" (199.4cm x 101.6cm)
Model 344	Torso Blanket 43" x 40" (109.2cm x 101.6cm)
Model 542	Sterile Cardiac Blanket 72" x 50" (183cm x 127cm)
Model 145	Warming Tube 28" x 66" (71.1cm x 167.6cm)