

MAXI-THERM®

SINGLE-PATIENT USE HYPER-HYPOTHEMRIA BLANKETS



The Maxi-Therm® blankets are durable for controlled warming and cooling for patients with longer stays. They have a random flow design which results in even temperature distribution. The design is so efficient the blanket can be filled with the patient in place, allowing for the blanket to be used on top or under the patient.

Maxi-Therm® Product Description

Model 276	Adult - 24" x 60" (61cm x 152.4cm)
Model 274	Pediatric - 22" x 30" (55.9cm x 76.2cm)
Model 273	Infant - 12" x 18" (30.5cm x 45.7cm)

**must purchase 286 Hose separately to use with equipment