

MICRO-TEMP® LT

LOCALIZED HEAT THERAPY SYSTEM

The Micro-Temp® LT consistently delivers a therapeutic level of localized heat using a pump, connecting hose and localized temperature therapy pad. A variety of safety features are incorporated into the Micro-Temp® LT. These features include audible and visual alarms for patient and product protection.

Additional protection includes an audible and visual alarm if temperature exceeds 1 °C over setpoint. The new Micro-Temp® LT will also shut down in the event the unit tips over or when the water level is too low.



Product Description

Model 749	100V, 115V & 220/240V (Not Available in Europe)
Dimensions	9.8" W x 8.2" D x 9.1" H (24.9cm x 20.8cm x 23.1cm)
Water Temp. Range	20° C - 42° C (68° F - 107.6° F) HEATING ONLY