

NORM-O-TEMP®

HYPERTHERMIA SYSTEM



The Norm-O-Temp® is a whole body system that allows conductive warming therapy to be administered. It can be used in the operating room, pre-op, recovery, or emergency department and when paired with the Gelli-Roll® makes the perfect Temperature Management device for your OR.

Product Description

Model 111Z	115v
Model 111W	230v
Dimensions	9" W x 15" D x 18.5" H (22.9cm x 38.1cm x 46.9cm)
Temperature Settings	20° C - 42° C (68° F - 107.6° F) HEATING ONLY